

Replacing Gearbox

- 1) Move the chair forward and upright for easy access
- 2) Turn off the power switch, then unplug the power cord
- 3) Lift up the pillow & cover cushion
- 4) unzip two(2) zippers to expose the massage mechanism
- 5) Remove back cover panel
- 6) Locate four(4) screws from gearbox, then remove 4 screws (Screws are behind the strap)
- 7) By using your finger, dislocate the dc belt from gearbox, then remove the gearbox from the back of the chair
- 8) Replace with new gearbox, then screw four(4) screws back to mount onto the channel frame
- 9) Install DC-belt to gearbox pulley
- 10) Close the back cover panel

